

Urgent® PC Treatment Schedule

Urgent® PC therapy begins with 12 treatments (30 min. each).

After your 12th treatment, you and your doctor will discuss whether you would like to continue Urgent PC therapy. For most patients, symptom improvements are sustained with a refresher treatment about once a month.

SESSION	APPOINTMENT
Week 1	
Week 2	
Week 3	
Week 4	
Week 5	
Week 6	
Week 7	
Week 8	
Week 9	
Week 10	
Week 11	
Week 12	
Follow up visit	(Schedule 3-4 weeks after 12 th treatment)
Follow up visit	
Follow up visit	
Follow up visit	
Follow up visit	
Follow up visit	
Follow up visit	



Non-drug, non-surgical
OAB Treatment
For men and women



www.AdvUro.com

How does the Urgent® PC PTNS Therapy work?

This therapy is a gentle and gradual way of modifying the OAB signals.



During treatment, mild electrical impulses:

1 Enter through an acupuncture-type needle placed near ankle



2 Travel up tibial nerve in the leg

3 Reach nerves responsible for bladder control



Why Urgent® PC PTNS Therapy

- Studies show that most patients can achieve long-lasting relief of their symptoms
- After Urgent® PC, most patients go to the bathroom less and have fewer accidents
- More than one million treatments since 2003
- Extremely low-risk and comfortable for most patients

Urgent® PC PTNS Therapy schedule

- 30 minute in-office treatments during which you can read, play games or catch up on email
- 12 weekly sessions to determine your ideal response. Most patients can expect to see some improvements by about 6 treatments but it can take up to 12 weeks to see results
- This is followed by periodic maintenance sessions for long-term relief. These can be done around your schedule so you can continue enjoying the freedom of your improved bladder control



Urgent PC

NEUROMODULATION SYSTEM

Treatment for Overactive Bladder

The Urgent PC SYSTEM

Urgent PC stimulator - delivers electrical impulses

Needle electrode - transfers impulses to the tibial nerve

Lead set - connects adhesive pad, needle electrode and stimulator



1 Get Comfortable

During treatment you will sit with your foot supported with a thin needle electrode inserted near your ankle. Impulses from the stimulator travel through the needle electrode, along your leg and to the nerves in your pelvis that control bladder function.



2 Needle Placement

Your clinician will press on your ankle to find the right place to insert the needle electrode and then clean that spot. You may feel a tap of the needle as it pierces your skin and is inserted.



3 Connect Simulator

Your clinician will attach a round patch to the bottom of your foot and a clip onto the needle electrode. These connect you to the stimulator.

